

Precision™XS

Automatic pipetting device for MICRONAUT plates

▼ What are the features of Precision™XS (a product of BioTek Instruments)?

MC Farland tubes are diluted automatically. The dispenser transfers this dilution with the necessary media into the microtiter plate.

▼ How fast works Precision™XS?

Pipettes up to 50 microtiter plates (up to 200 tests) within 2 hours.

▼ Further advantages of Precision™XS?

No need of additional Bouillon tubes or MICRONAUT reservoirs.

▼ Combined with BioStack™ the automated liquid handling throughput of 50 microtiter plates can be achieved.

▼ Technical data

Measure: 51 x 76 x 56 cm
Weight: 29 kg
Software: Precision Power™ PC software
Dispenser: 1-channel- and 8-channel-pipettes, 8-channel dispenser

